

## SarcomaFusion test




The **Genexpath SarcomaFusion** solution can highlight 140 fusion transcripts responsible for sarcomas. It is applicable to low-quality RNA samples extracted from sections of FFPE biopsies, fixed and included in paraffin, obtained in the clinic. A needle biopsy is also possible to obtain a sample of sufficient quality. Data are generated using a next-generation sequencer and require only  $10^5$  reads per sample.

They are then processed using our RT-MIS software which provides details of any transcripts detected and the associated bibliography.

### Using the test

The **GENEXPATH SarcomaFusion** test is based on a ligation-dependent RT-PCR (LD-RT-PCR) method. This semi-quantitative technique makes it possible to simultaneously evaluate the expression levels of a large number of genetic markers and more precisely for this test, chromosomal translocations using pairs of oligonucleotide probes specific to each of these markers.



### A simple and fast protocol

From a total RNA extract, four steps are sufficient to obtain the libraries.

1. A reverse transcription step.
2. A step of hybridization of oligo-nucleotide probes.
3. A ligation step.
4. A PCR amplification step.
5. Then sequencing the libraries

No purification is necessary until the libraries are obtained, which limits the loss of material and ensures a very good sensitivity to this technique. In addition, the genetic sequences targeted by the probes are particularly short (between 40 and 60 bases) which ensures a very good robustness against the degradation of RNAs.

LD-RT-PCR is therefore a particularly suitable approach for the analysis of difficult biological samples such as tissue biopsies fixed and included in paraffin.

For each sample, about  $10^5$  sequences are sufficient to obtain an analyzable expression profile, which makes it possible to test a large number of samples in parallel on the same sequencing FlowCell.

To optimize costs, **GENEXPATH SarcomaFusion** libraries can also be loaded at the same time as other sequencing libraries, generated by other methods.

### **Post-PCR analysis based on dedicated software**

Once the sequencing is complete, the FASTQ file can be uploaded to the RT-MIS platform which, after a few minutes of analysis, delivers a file containing the detected merge transcripts.

RT-MIS also provides a bibliography for detected transcripts and thus offers a complete solution for cancer researchers.

<b>Handling duration</b>	5hrs30'
<b>Actual working time</b>	1hr-1hr30'
<b>Type of nucleic acid</b>	RNA
<b>Input quantity</b>	Between 50 and 500ng of RNA in a volume of 2µl
<b>Type of cancer</b>	Sarcoma
<b>Contents of the reagent kit</b>	Probes targeting more than 140 fusion transcripts, barcodes, sequence primers
<b>Method</b>	Ligation dependent RT-PCR
<b>Description</b>	Detects more than 140 transcripts of sarcoma-associated fusion in 1 single analysis.
<b>Equipment compatibility</b>	MiSeq, NextSeq 500, NextSeq 550 Illumina®
<b>Type of samples</b>	Tissue biopsies fixed and included in paraffin
<b>Technology</b>	Next Generation Sequencing



*In vitro diagnostic medical device according to Directive (EU) 98/79/EC*



*For in vitro diagnosis use only. Please read the user manual before use.*

## Contact

Genexpath, 113 avenue des martyrs de la résistance 76100 Rouen, France

Phone : +33 (0)2 78 08 98 69

E-Mail : [contact@genexpath.com](mailto:contact@genexpath.com)

Website : [www.genexpath.com](http://www.genexpath.com)