

In vitro protocol – in brief test SarcomaFusion

1. Reverse Transcription

Reagents:

- 5X Vilo reaction mix
- 10X super script
- Nuclease free water
- RNA samples

Steps :

- Traw reagents
- Prepare reverse transcription mix
 - 5X Vilo reaction mix (1 μ L)
 - Nuclease free water (1 μ L)
 - 10X super script (0,5 μ L)
- Deliver 2,5 μ L of mix per tube
- Add 2,5 μ L of RNA sample
- Vortex and centrifuge
- Put samples in thermocycler and start program 1 step 1
- Once the program is over and the block cool down at 4°C, put out tubes
- Centrifuge
- Put tubes on ice or cooling block



Then proceed directly to step 2 or store cDNA products between -30°C and -15°C.

2. Hybridization of probes

Reagents:

- SALSA MLPA Buffer
- GEP- SFPM** Probe Mix

Steps :

- Traw reagents
- Prepare hybridization mix
 - Salsa MLPA Buffer (1,5 μ L)
 - GEP- SFPM** Probe Mix (1,5 μ L)
- Vortex and centrifuge
- Deliver 3 μ L of mix per cDNA tubes
- Centrifuge
- Put samples in thermocycler and start program 1 step 2

3. Ligation

Reagents:

- SALSA Ligase Buffer A
- SALSA Ligase Buffer B
- SALSA Ligase 65
- Nuclease free water

Steps :

- Traw reagents
- Prepare ligation mix
 - Nuclease free water (25 μ L)
 - Salsa Ligase Buffer A (3 μ L)
 - Salsa Ligase Buffer B (3 μ L)
- Vortex and centrifuge
 - Salsa Ligase 65 (1 μ L)
- After 1h of hybridization, deliver 32 μ L of mix per cDNA tubes
- Start program 1 step 3



Then proceed directly to step 4 or store the ligation products between -30°C and -15°C.



After this step, do not store the products at high temperatures (e.g. 4°C or at room temperature) in order to avoid non-specific ligations which could result from a residual activity of the enzyme.

4. Amplification and incorporation of barcodes and adaptators

Reagents:

- GEP-BC-xxx Barcodes
- GEP-BCC-xxx GAPDH Barcodes
- Red'y' Star PCR Mix
- Nuclease free water

Steps :

- Traw reagents
- Prepare PCR mix
 - Red'y' Star PCR Mix (12,5 µL)
 - Nuclease free water (5,5 µL)
- Vortex and centrifuge
- Deliver 18 µL of mix per wells of PCR plate
- Add 5 µL of ligation product in each well
- Add 2 µL of barcode (GEP-BC-xxx or GEP-BCC-xxx following the test) in each well



Use different barcodes BEP-BC-xxx for each tested samples.

- Put samples in thermocycler and start program 2



Then proceed directly to step 5 or store the PCR products between -30°C and -15°C.



Do not store these products for prolonged periods at high temperatures (e.g. 4°C in the thermocycler or at room temperature).

5. Purification and dosage of sequencing libraries

Reagents:

- 100% Ethanol
- Nuclease free water
- AMPure XP Magnetic beads
- TE Buffer
- Qubit dsDNA HS Assay

Steps :



Ensure beads are completely re-suspended before use.

- Purify 25 µL of PCR products with 45 µL AMPure XP
- Add 50 µL of TE buffer to each well of the reaction plate



After purification, libraries can be stored between -30°C and -15°C before sequencing.

- Dose 10 µL of each library by fluorimetry



This rapid protocol completes the notice. It does not dispense with the complete reading of the notice.