

## Sequencing protocol – in brief

### test SarcomaFusion

### Illumina MiSeq System

#### Reagents:

- Illumina sequencing reagents

#### 1. Dilutions and denaturation

##### Steps :

- Dilute each library at 2-4 nM concentration
- Pool libraries in equivolume  
If other libraries are sequenced, adjust pool concentrations and then combine them
- Denature and dilute final pool at 8-10 pM loading concentration

#### Reagents:

- GEP-SP-001** primer
- GEP-SP-002** primer
- Illumina sequencing reagents

#### 2. Preparation of sequence primers

##### Steps :

If only pool of SarcomaFusion libraries

- Dilute primers
  - GEP-SP-001** (3  $\mu$ L)
  - GEP-SP-002** (3  $\mu$ L)
  - HT1 Buffer (594  $\mu$ L)

- Load 600  $\mu$ L into reservoir #18 of the flowcell

If pool of SarcomaFusion libraries is associated with others libraries

- Pipet 600  $\mu$ L from reservoir #12
- Add
  - GEP-SP-001** (3  $\mu$ L)
  - GEP-SP-002** (3  $\mu$ L)
- Load all volume into reservoir #18 of the flowcell

#### 3. Preparation of the injection sheet

- If the GENEXPATH SarcomaFusion library is sequenced alone, create the injection sheet to generate the FASTQs, providing for 120 cycles in read 1.
- If the GENEXPATH SarcomaFusion libraries are combined with other sequencing libraries, generate the injection sheet using the usual parameters, without entering the GENEXPATH SarcomaFusion samples.
- Specify the use of custom during run setup (With Local Run Manager, on the Create Run page. In manual run mode, on the Run Setup screen).



**In all cases, ensure that reading in read 1 is done with a minimum of 120 cycles and that the Custom Primer for Read 1 box is selected.**

#### 4. Sequencing start

Initiate sequencing following the procedure described in the MiSeq System Illumina Guide.



**This rapid protocol completes the notice. It does not dispense with the complete reading of the notice.**